



# Large Particle Flow Cytometry Platform for High-Throughput Analysis, Brightfield Imaging, and Precision Dispensing of Intact HT-29 and HT-1080 Spheroids

Megan S. Harrison<sup>1</sup>, Deborah R. Frenkel<sup>1</sup>, Yongwoon Kim<sup>1</sup>, Rock Pulak<sup>1</sup>, Hilary Sherman<sup>2</sup>, Ann Ferrie<sup>2</sup>  
 1 Union Biometrica, Inc., Holliston, MA USA, 2 Corning Life Sciences, Kennebunk, ME USA

## Introduction



**Union Biometrica, Inc.** develops instrumentation that extends the principles of flow cytometry to large multicellular structures and whole organisms.

These platforms enable researchers to quantitatively analyze and selectively isolate intact biological systems including:



- 3D cell culture models (organoids, spheroids)
- Small model organisms (*C. elegans*, zebrafish)

**Spheroids:** multicellular 3D aggregates that recapitulate key aspects of tissue architecture and gradients  
**Organoids:** stem cell-derived systems that self-organize into structures resembling native tissue

**Why 3D models matter:**


- More physiologically relevant than traditional 2D culture
- Increasingly used in drug discovery, disease modeling, and regenerative medicine
- 3D models are central to the shift towards human-relevant, non-animal models and new approach methodologies (NAMs)
- FDA draft guidance mandates that NAMs—including organoids and other 3D models—must meet four criteria: fit-for-purpose, biological relevance, technical characterization, and reproducibility

## Materials and Methods



GFP-positive HT-29 and non-fluorescent HT-1080 cell lines were cultivated as roughly uniform, ~200-micron spheroids using Corning® Elplasia® 12K flasks. HT-29 spheroids and HT-1080 spheroids were mixed 1:1, then data were acquired from mixtures of heterogeneous, intact spheroids using the COPAS Vision equipped with a 1000-micron flow cell. HT-29 spheroids were distinguished from HT-1080 spheroids based on Green Fluorescence. GFP-positive HT-29 spheroids of uniform size were automatically dispensed singly into the wells of a 96-well plate. Brightfield imaging, Profiler, and Forward Scatter provide additional means to distinguish HT-29 from HT-1080 spheroids.


## Flow Cytometry of Intact Spheroids



Measurements include:


- Size (TOF - Time Of Flight)**
- Optical Density (Extinction)**
- Fluorescence**
- Scatter**

## Population-level Batch Analysis Data




- HT-29 spheroids and HT-1080 spheroids make up 51.1% and 45.9% of all spheroids respectively
- HT-29 spheroids have a mean diameter of **230 microns**
- HT-1080 spheroids have a mean diameter of **211.8 microns**

## Precision Plate Dispensing of Intact Spheroids




- Dispensed one GFP+ HT-29 spheroid per well of a 96-well plate
- Filled plate in **~2 minutes** with **100% accuracy**

## Brightfield Imaging in-flow




- HT-29 and HT-1080 spheroids can be distinguished based on the appearance of brightfield images
- HT-29 spheroids appear dark, compact, and smooth
- HT-1080 spheroids appear light, loose, and rough
- Image analysis provides information on spheroid size in microns or microns<sup>2</sup>

## Profiler™: Beyond the Integral Signal



- Profiler displays values for each measured parameter as they vary
- The Peak Height, Peak Width, and Integral signal (area under the curve) are captured for each parameter
- The Green Fluorescence profile for HT-29 spheroids is consistent with uniform localization of signal throughout the spheroid

## Label-free Parameters Distinguish Spheroids



Researchers may avoid fluorescent labeling in cases where viability or biological activity may be impacted

- HT-29 spheroids can be distinguished from HT-1080 spheroids based on differences in Forward Scatter

## Summary

- Large particle flow cytometry enables high-throughput analysis of intact 3D cellular models without dissociation.
- Combined cytometric measurements, Profiler data, and in-flow brightfield imaging provide orthogonal characterization of spheroid populations in a single workflow
- Population-level analysis reveals heterogeneity in size, morphology, fluorescence, and optical density across intact spheroid models.
- Precision dispensing enables rapid isolation of selected intact spheroids into multi-well plates for downstream assays and workflow automation.
- These capabilities support scalable, information-rich analysis and handling of complex 3D biological structures.

